

勝特力材料 886-3-5753170
勝特力电子(上海) 86-21-34970699
勝特力电子(深圳) 86-755-83298787
Http://www.100y.com.tw

HAMAMATSU
PHOTON IS OUR BUSINESS

STABILIZED LIGHT SOURCES FOR PHOTOMULTIPLIER TUBES



L 11494




L 11416


L11494 Series
Image during use

- ① L11494
- ② Photodetector

TACCC0150EB



L11416 Series
Image during use



TACCC0149EA

SPECIFICATIONS

L11416 Series

(at +25 °C)

Parameter	L11416-470	L11416-525	L11416-555	L11416-590	Unit
Peak Emission Wavelength	465	522	555	592	nm
Spectral Half Width	26	35	30	18	nm
Light Emission Power	1 ± 0.5 *1				pW
Emission Stability (0 °C to +50 °C)	±2 Max.				%
Battery	Button battery SR41 or equivalent				—
Battery Service Life	Min.	24			h
Operating Ambient Temperature / Humidity *2	0 °C to +50 °C / below 85 %				—
Storage Temperature / Humidity *2	-20 °C to +60 °C / below 85 %				—
Weight	13				g

*1: Adjusted by using a 22 mm effective diameter photon counting head placed at a position 5 mm away from the L11416.

*2: No condensation

L11494 Series

(at +25 °C)

Parameter	L11494-430	L11494-470	L11494-525	L11494-660	Unit
Peak Emission Wavelength	428	465	522	660	nm
Spectral Half Width	65	26	35	30	nm
Light Emission Power	High: 1 ± 0.5, Low: 0.01 ± 0.005 *3				pW
Emission Area	φ7.0				mm
Emission Stability (0 °C to +50 °C)	±2 Max.				%
Battery	Button battery SR41 or equivalent				—
Battery Service Life	Min.	10	24		h
Operating Ambient Temperature / Humidity *4	0 °C to +50 °C / below 85 %				—
Storage Temperature / Humidity *4	-20 °C to +60 °C / below 85 %				—
Weight	23				g

*3: Adjusted by using an 8 mm effective diameter photon counting head placed in direct contact with the L11494.


*4: No condensation

OVERVIEW

The L11416 and L11494 series are LED light sources mainly designed for adjusting the sensitivity of PMT (photomultiplier tubes) and PMT modules. The L11416 series is a test tube type easily mounted in lab devices such as for blood sampling, while the L11494 series is a plate type designed to be placed on a sample stage. Stable continuous light of approximately 1 pW is always output because the LED emission power is controlled by monitoring with a photodiode.

INTERNAL DIAGRAM


STABILIZED CONTINUOUS LIGHT




CHARACTERISTICS

EMISSION SPECTRUM CHARACTERISTICS


●L11494-430




●L11416-470/L11494-470




●L11416-525/L11494-525




●L11416-555




●L11416-590




●L11494-660



BATTERY SERVICE LIFE CHARACTERISTICS (Typical example)




TEMPERATURE STABILITY




DIMENSIONAL OUTLINES (Unit: mm)

L11416 Series



TACCA0310EA

L11494 Series



TACCA0311EA

勝特力材料 886-3-5753170
 勝特力电子(上海) 86-21-34970699
 勝特力电子(深圳) 86-755-83298787
[Http://www.100y.com.tw](http://www.100y.com.tw)

Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office. Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2012 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: Hamamatsu Corporation, 360 Foothill Road, P. O. Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road Welwyn Garden City Hertfordshire AL7 1BW, United Kingdom, Telephone: 44-(0)1707-294888, Fax: 44(0)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Thorshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia: S.R.L.: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Road North, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

TACC1057E01
DEC. 2012 IP