

Ring Fluorescent Lamp for Stereo Microscope

Our Fluorescent Lamp Illuminator offers shadow-free and non-Directional illumination by fitting lamp onto microscope objective lenses. It can be used with all kinds of Microscopes. It is ideal lighting for long working without eye fatigue.

Color temperatures, RED / GREEN / BLUE / Yellow filter glass made, YEX / Blacklight are available for special requirements.

?Application (for industrial purpose)

IC Inspection, Jewel Inspection, and Production control by ITV,

(for science)

physiological / pathological inspection

agricultural research




勝特力材料 886-3-5753170
 勝特力电子(上海) 86-21-54151736
 勝特力电子(深圳) 86-755-83298787
[Http://www.100y.com.tw](http://www.100y.com.tw)

The type of standard lamp and its characteristics

Lamp Type No.	Color Temperature K	Spectrum Distribution A	Peak Spectrum A	Application
80B-WW 92B-WW WarmWhite	3200 3950	4000 7200	5800	Standard White for usual application
80B-DL 92B-DL DayLight	5300 8000	4000 7200	4600	Bluish color, fit to be used for inspection of metal surface

Warm White (WW)



Daylight (DL)

